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**Warning notices:** The product must be cleaned, disinfected and sterilised before first use as well as after any further use, according to the enclosed instructions

GB

**Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004**

Hersteller: HAGER & WERKEN GmbH & Co. KG  
Ackerstraße 1, 47269 Duisburg, Germany  
T +49 (203) 99269-0, F +49 (203) 299283  
www.hagerwerken.de

Product:  
Fit Clamp Forceps (REF 355 333)  
Fit Punch (REF 355 321)

**Limitation of reprocessing:** Frequent reprocessing has only got a small impact on this instrument. The instrument's shelf life is mainly determined by wear and damage by use.

**Place of handling:** Clean surface with a disposable cloth or paper.

**Storage and transport:** No special requirements.

**Preparations for cleaning:** Disassemble product into its components.

**Cleaning:** Only use suitable cleansers, carefully following the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C in mechanical cleaning.

**Disinfection:** Only use suitable disinfecting solutions, carefully following the instructions of the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C in mechanical disinfection.

**Maintenance:** No special requirements.

**Packaging:** Standardised packaging material may and should be used for sterilisation processes.

**Sterilisation:** Vapour sterilisation at 134 °C for 5 minutes, vapour sterilisation at 121 °C for 15 minutes.

**Control/Functional check:** Sight check on damages, wear, deformation.

**Storage:** No special requirements.

## MEDICAL DEVICE

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine controls are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

2015081701 Rev6 (2020-10)

HAGER & WERKEN GmbH & Co. KG · Ackerstraße 1 · 47269 Duisburg, Germany  
T +49 (203) 99269-0 · F +49 (203) 299283 · www.hagerwerken.de