MIRASUC® ADAPTER II

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GB

Processing instructions according to EN ISO 17664:2021

Product: Mirasuc Adapter II, white, 6/11mm

REF 254 132

Manufacturer: HAGER & WERKEN GmbH & Co. KG, Ackerstraße 1, 47269 Duisburg, Germany

Warnings: Warning! The product must be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every subsequent use. Disinfection alone is not enough. Thorough cleaning and disinfection are essential prerequisites for effective sterilization. The processing should be started as early as possible, at no later than 2 hours after use. During the processing, the product should not be unnecessarily exposed to moisture or humidity. Please also observe the legal regulations applicable in your country as well as the hy-giene regulations of your medical practice and/or the hospital. **Limitations in the processing:** Frequent reprocessing has only minor effects on this product. The end of the product's useful life is primarily determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced at an early stage. The user is solely responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product. **Personal protective equipment:** For reasons of occupational safety and to minimize infections/cross-infections, suitable personal protective equipment (protective clothes, protective gloves, protective goggles and nose/mouth protection) must be used during the entire processing procedure.

INSTRUCTIONS:

INITIAL TREATMENT AT THE PLACE OF USE

Remove surface contamination using a disposable cloth/paper towel while still at the place of use. Rinse the products with water (at least drinking water quality) no later than 2 hours after use. Do not allow residues or contaminations of any kind to dry on the product. For the initial treatment at the place of use, do not use any aldehyde-containing or alcohol-containing agents as they may cause protein fixation.

PREPARATION BEFORE CLEANING

If possible, disassemble the product into its individual components. Inspect for damage and wear.

CLEANING & DISINFECTION: BY MACHINE (RECOMMENDED)

Only perfectly maintained, checked, calibrated and released cleaning and disinfection equipment according to EN ISO 15883 with valid test mark (CE mark) may be used. Make sure that during machine cleaning and disinfection of several different products/individual components, the cleaning device is not overfilled and that the products/individual components are not in contact with each other. We recommend completing the cleaning and disinfection cycle according to the specifications of the manufacturer of the device and the cleaning agent and/or disinfectant. Make sure that the cleaning agent is compatible with the products and disinfectants that might be used. The pH of the cleaning agent and disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Application with the following agents is recommended: Dr. Weigert MediClean forte 0,5 %. The instructions for use of the device manufacturer and the manufacturer of cleaning agent and disinfectant must be strictly observed. In particular, the concentrations and contact times have to be observed.

Description of the recommended cleaning process: Pre-cleaning: The pre-cleaning should be completed in a water bath (at least drinking water quality) without other cleaning additions using a cleaning brush below the water surface to achieve both sufficient cleaning of the products/individual components while avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; it should, however, not be less than 1 minute.

Cleaning tools for pre-cleaning: Mirasuc brush REF 605 216

Information about the cleaning and disinfection device: Type: Miele G7881

Description of the recommended cleaning and disinfection: Means: Please use the following flushing nozzles as adapters and before the cleaning and disinfection process; attach the products onto the flushing nozzles. For the exact application of the flushing nozzles, please refer to the instructions of the relevant manufacturer. Flushing nozzles A 833, Miele & Cie KG, Gütersloh, Flushing nozzles A 816, Miele & Cie KG, Gütersloh. **Cleaning program:** VarioTD

Cleaning temperature: 50 - 56.43 °C **Hold time:** 10:50 min **Disinfection temperature:** 65 - 94.4 °C **Hold time:** 09:52 min **Rinsing technique:** Fresh water system/demineralized water **Cleaning agent:** Dr. Weigert MediClean forte 0,5 % **Neutralizer:** Dr. Weigert MediClean forte 0,5 %

The water should be of drinking water quality. Flushing and drying should be completed by machine. The machine cleaning should not exceed a temperature of 60 °C. Ensure that visible contamination has been completely removed. If there is still visible contamination after the cleaning and disinfection process, the cleaning & disinfection must be repeated. Before further processing, the product must be free from any residues and dry.

CLEANING: MANUAL

Manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both sufficient cleaning of the products/individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should

be based on the level of contamination of the product/individual component; it should, however, not be less than 1 minute. Afterwards, rinse the products/individual components with water (at least drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process: **Pre-cleaning:** Remove surface contamination using a disposable cloth **Cleaning:** Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution; remove contamination by means of cleaning brushes. **Cleaning agent:** ID 213, company Dürr Dental, concentration: 2 % **Cleaning tools:** Mirasuc brush, REF 605 216 **Rinsing:** Rinse with water (at least drinking water quality) for 20 sec. **Drying:** Dry at room temperature.

Manual cleaning should not exceed a temperature of 45 °C. Ensure that visible contamination has been completely removed. If there is still visible contamination after cleaning, repeat the cleaning process.

DISINFECTION: MANUAL

For manual disinfection, please use only approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products/individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not allow them to be in contact with each other. The instructions for use of the disinfectant manufacturer must be strictly observed. In particular, the concentrations and contact times must be observed. Application with the following agents is recommended: Dürr Dental ID 213 instrument disinfection.

Description of the recommended disinfection process: Completion of proper cleaning preparation and manual cleaning **Disinfection:** Insert the instrument into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution **Disinfectant:** ID 213, company Dürr Dental, concentration: 2 % **Cleaning tools:** Mirasuc brush, REF 605 216 **Contact time:** 5 min **Rinsing:** Rinse with water (at least drinking water quality) for 20 sec. **Drying:** Dry at room temperature.

Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Please observe thermal restrictions that might apply considering the disinfectant manufacturer's specifications. Afterwards, rinse the products/individual components with water (at least drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry.

DRYING

If you process the products manually, the products/individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature. When machine processing, drying should be completed directly in the devices using hot air. With machine drying, the temperature of 100 °C should not be exceeded.

MAINTENANCE, CONTROL AND CHECK

The product does not have special requirements for maintenance. Re-assemble disassembled products/individual components. You must always carry out a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product/individual component is not visibly clean, the entire processing must be repeated or the product/individual component must be properly disposed of.

PACKAGING

The product is to be packed in suitable and standardized transparent sterile packaging (sterilization bags) and sealed. Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products/individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION

Only tested steam vacuum autoclaves may be used. Make sure that during the sterilization of several different products/individual components, the autoclave is not overfilled and that the products/individual components do not come into contact with each other. The following sterilization cycles may be completed: Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 15 minutes. To dry the products/individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer's instructions for use.

STORAGE

To maintain the sterility, the products should be stored in standardized sterilization bags in a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed once again before their use.

ADDITIONAL INFORMATION

Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling. The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded. The manufacturer and the competent authority of your member state shall be immediately notified all serious incidents occurring in connection with the product. For the disposal of the products, there are regionally differing recommendations and regulations. Ask your competent disposal company for the current regulations in your region.

CONTACT TO THE MANUFACTURER

Phone: (+49) 203/99269-0 Fax: (+49) 203/299283 Email: info@hagerwerken.de

The instructions specified above are validated by the medical device manufacturer as suitable for the preparation of a medical device for its reuse. The processor is responsible for ensuring that the actually completed processing with the equipment used, materials and personnel in the preparation facility achieves the desired result. To this end, verification and/or validation and routine monitoring of the procedure are required.

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