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Warning notices: The product must be cleaned, disinfected and sterilised before first use as well as after any further use, according to the enclosed instructions.

GB

Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: HAGER & WERKEN GmbH & Co. KG
Ackerstraße 1, 47269 Duisburg, Germany
T +49 (203) 99269-0, F +49 (203) 299283
www.hagerwerken.de

Product: Emmenix Filmholder
REF 656 300
valid as of 07/2018 (Rev: 3)

Limitation of reprocessing: max. 20 x

Place of Handling: Clean surface with a disposable cloth or paper towel.

Storage and Transport: Storage and Transport.

Preparations for Cleaning: n/a

Cleaning: Only use suitable cleansers according to the respective manufacturer instructions. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C for mechanical cleaning.

Disinfection: Use suitable disinfecting solution according to the respective manufacturer instructions. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C for mechanical disinfection.

Maintenance: No special requirements.

Packaging: Standardized packing material for sterilisation can be used (recommendation).

Sterilization: Vapor sterilization at 134 °C for 5 minutes or vapor sterilization at 121 °C for 15 minutes.

Control/Functional Check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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Fabrikant: HAGER & WERKEN GmbH & Co. KG