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REF 452 403 bis REF 452 411, REF 452 415, REF 452 416

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Instructions for Use

REF 452 403 to REF 452 411, REF 452 415, REF 452 416

Intended use: Electrosurgery accessories must only be used by qualified healthcare professionals! This instruction does not replace the reading of the instructions for use of the electrosurgery device and further accessories used. HF electrodes are intended for open surgery and are used for cutting and coagulating biological tissue. They are not intended for having direct contact with the heart or the central circulatory system. HF electrodes are connected to the monopolar output of the electrosurgical generator via an HF cable or an HF handle. The electrodes are activated by means of a foot switch or an HF handle.

Indication: Monopolar electrode for cutting and coagulating soft tissue during surgical-dental interventions, for example:

- cutting the frenulum of the lip
- opening acute abscesses
- in case of parodontal diseases
- for removing tumours and fistulas
- for exposing retained teeth

Contraindication:

- Restrictions for use in connection with cardiac pacemakers and in-vivo cardiac defibrillators
- Flammable vapours and liquids – explosion hazard
- Do not use in immediate proximity to bone or the root of a tooth
- Avoid contact with metal

Undesired side effects:

- The contact with metallic objects results in undesired remote effects

Safety advice:

Caution: The max. peak voltage of 1000 VP or 1 KvP must not be exceeded!

Before first use, new HF electrodes must be cleaned and sterilised according to a validated procedure (DIN EN ISO 17665). A functional test must be carried out before every use. It must be made sure that the HF electrode is securely inserted into the handle. This must be performed with care in order to avoid damage to the HF electrode and/or injuries of the patient or the surgical staff. The tip of the electrode may be damaged when excess force is applied. The HF electrode must not be activated while it is in contact with metallic objects and/or optics. You must make sure that no flammable material is in immediate proximity during the complete use since otherwise there will be an explosion hazard. After the electrosurgical current has been switched off, the tip of the electrode may still be hot and cause burns. The use of HF current may cause damage to cardiac pacemakers and in-vivo cardiac defibrillators, therefore persons affected must seek advice from a cardiologist before the intervention. Make sure to insulate metallic restorations using plastic matrices as well as to use dental mirrors, aspirators, cheek and lip retractors made of non-conducting materials. When exposing cavity and crown borders, make sure not to go below a biological width of 3 mm.

Performance data: The HF generator should provide the following output power:

- cutting mode 50 watts at 600 Ω
- coagulation mode 30 watts at 1200 Ω

Warning: Tips of electrodes may cause injuries!

- After use, tips of electrodes may be very hot and cause burns!
- Never place electrodes on the patient or in their immediate proximity!
- Do not use in the presence of flammable or explosive substances!
- If the electrode is used together with rinsing instruments, if possible, use a non-conducting rinsing fluid!
- Deep incisions must not be made in one cut, adhere to the corresponding cool-down time in order to avoid over-heating or burning the deep tissue resulting in wound healing disorders.

Recommended devices/Exception of Siemens devices: We recommend the use of generator units manufactured by Gebrüder Martin, Berchtold, Erbe und Meyer-Haake. Compatibility with the handle for 1.6 mm shafts is required. To connect the electrosurgical accessories (e.g. handle, electrodes), the correspondingly compatible cables according to the manufacturer's instructions must be used.

Warning: Our electrodes cannot be used together with devices/units manufactured by Siemens, in particular with the M1 model since the isolation becomes very hot and may cause burns.

Before use: Before every use, the insulation of the HF electrodes must be checked for indentures or damage. HF electrodes that present damage or indentures must not be used: the electrode must be checked for intact insulation, cleanliness and integrity. Only use flawless and sterilised products! It must be made sure that there is a safe and perfect connection between electrodes, handle and generator before start of operation. Therefore, insert the electrode connection into the dedicated electrosurgery handle completely and according to instructions, no metal parts of the electrode shaft may be visible any more. The maximum permissible operating voltage of the handle must not be exceeded.

Warning: Connect the electrode to the electrosurgery device only when it is switched off or in standby mode. Failure to observe this instruction may result in burns or electric shock! Make sure to place the neutral electrode correctly!

During use: Always work using the lowest performance setting for the desired surgical effect.

Cleaning and sterilisation: Not sterile, clean and sterilise before initial and all further use.

Notes: Municipal water to be used must comply with the Council Directive 98/83/EC of 03 November 1998 on the quality of water for human consumption.

This instruction for reprocessing states the cleaning agents and disinfectants used for validation. The person who performs the reprocessing must assume all responsibility if they use an alternative cleaning agent and disinfectant (RKI- or VAH-listed). Due to the product design, the materials used and the intended use, no defined limit of reprocessing cycles that can be performed can be established, however, corresponding tests confirm that there is no functional impairment of the electrodes after 20 reprocessing cycles. The service life of the instruments is determined by their function and their gentle use. Naturally, instruments for electrosurgery are subject to increased wear and tear depending on the type and duration of use.

Instruction for reprocessing

Preparation for cleaning: Remove the HF electrodes from their packaging. They must be placed in a container/device intended for cleaning/sterilisation.

Preparation and transport: Remove coarse dirt from the instruments directly after use.

Information on cleaning in the CSSD: Immediately after use, immerse the products in cold water (< 40 °C). Do not use any fixing agents or hot water (> 40 °C) since this results in a fixation of residues and may impair the cleaning result. The electrodes must be stored and transported to the place of reprocessing in a closed container in order to avoid damage to the sensitive electrodes and environmental contamination.

Manual precleaning: Rinse the product using cold municipal water (< 40 °C) until all visible contamination has been removed. Use a soft brush to remove any tenacious dirt. Use a high-pressure water sprayer (or the like) to rinse cavities, lumens thoroughly (> 30 sec.) with cold municipal water (< 40 °C).

Manual reprocessing

Cleaning

1. Immerse the products in an alkaline detergent (e. g. 0.5 % neodisher® LM2) in an ultrasonic bath for a time of exposure to ultrasonic waves of 10 min. and a frequency of 35 kHz. For doing so, follow the instructions provided by the manufacturer of the detergent.
2. Post-clean the products completely using a soft brush. Rinse any existing cavities and lumens thoroughly (> 30 sec.) using a high-pressure water sprayer (or the like).
3. Rinse the products under running municipal water to remove the detergent (> 15 sec.).

Disinfection

1. Immerse the products in an RKI- or VAH-listed disinfectant. For doing so, follow the instructions provided by the manufacturer of the disinfectant. It must be made sure that the disinfectant reaches absolutely all areas of the product.
2. The process is validated for the following disinfectant: 3 % Korsolex® Plus, 15 min.
3. Rinse the products (complete washing inside, outside and cavities) in purified water > 15 sec.

Drying: Manual drying using a lint-free disposable wipe. In order to largely avoid water residues in cavities, we recommend air blowing using sterile, oil-free compressed air. The product must never be heated over 140 °C.

Machine reprocessing (Miele G7835 CD)

Cleaning: Place the instruments onto a sieve tray on the carriage and start the cleaning process.

1. 1 min. precleaning with cold municipal water < 40 °C
2. Draining
3. 3 min. precleaning with cold municipal water < 40 °C
4. Draining
5. 5 min. cleaning at 55 °C ± 5 °C with 0.5 % alkaline detergent (0.5 % neodisher® Mediclean Dental)
6. Draining
7. 3 min. neutralisation (neodisher® Z, 0.1 %) with cold municipal water (< 40 °C)
8. Draining
9. 2 min. rinsing with purified water
10. Draining

The special instructions provided by the manufacturer of the cleaning machine must be followed

Disinfection: Machine thermal disinfection in cleaning and disinfecting device, meeting the national requirements as to the A₀-3000 value; > 5 min. at 92 °C ± 2 °C

Drying: Automatic drying according to the automatic drying procedure of the cleaning and disinfecting device 30 min. at 60 °C ± 5 °C. If required, subsequent manual drying using a lint-free cloth and air blowing of lumens using sterile oil-free compressed air.

Functional test: Before every use, the insulation of the HF electrodes must be checked for indentures or damage. HF electrodes that present damage or indentures must not be used: the electrode must be checked for intact insulation, cleanliness and integrity. If required, repeat the reprocessing procedure until the instrument is visually clean.

Sterilisation: Sterilisation of the products using a fractionated pre-vacuum procedure (in compliance with ISO 17665-1) respecting the respective national requirements. Sterilisation must be carried out using a fractionated pre-vacuum procedure with the following parameters.

- 3 pre-vacuum periods with a pressure of at least 3 bars +/- 0.5
- Sterilisation temperature of at least 134 °C
- dwell time: ≥ 5 min.
- drying time: at least 20 min.

Warning: Do not sterilise in hot air. Differing process parameters such as e.g. higher temperatures > 134 °C and a longer sterilisation time > 5 min. will reduce the service life of the electrode.

Information on the validation of reprocessing: The following test instructions, materials and machines were used for validation:

Cleaning agents:

- Alkaline detergent for manual reprocessing neodisher® LM2
- Alkaline detergent for machine reprocessing neodisher® mediclean Dental

Cleaning and disinfecting device: Miele G7835 CD

Cleaning programme: DES-VAR-TD

Neutralising agent: neodisher® Z

Steam steriliser: Tuttnauer EHS3870

See report for details:

- Manual reprocessing / Report of Zwisler Laboratory # 1612.1459
- Machine reprocessing / Report of Zwisler Laboratory # 1612.1455
- Sterilisation / Report of Zwisler Laboratory # 1612.2119

If the chemicals and machines stated above are not available, it is incumbent upon the user to validate their procedure accordingly.

Repair and modification

Defective HF electrodes must not be repaired. They must be replaced by new HF electrodes. Unauthorised modifications and repair measures are strictly forbidden and will result in a loss of the manufacturer's warranty. Any change in the product or failure to comply with these instructions for use result in the exclusion of liability by Otto Leibinger GmbH.

Packaging, storage, transport, handling

The HF electrodes must be stored in a clean and dry environment. They should be stored in a protective container with individual compartments or in special electrode containers. HF electrodes must be transported,

cleaned, maintained, sterilised and stored with utmost care. This applies in particular to fine tips and other sensitive areas. The operator must guarantee the maintenance of the sterile condition after the sterilisation process.

Returns

Returns will only be accepted if they are marked as “no hygienic risk” or “not decontaminated” and packed safely for transport.

Disposal

The HF electrodes, the packaging material as well as the accessories must be disposed of according to the applicable national provisions and laws.

Subject to change without notice.

Symbols

	Manufacturer
	CE label including the number of the notified body, mdc medical device certification GmbH, Kriegerstrase 6, 70191 Stuttgart/Germany
	Follow the instructions for use
	Production year/Batch number
	Not sterile
	Caution!

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HAGER & WERKEN GmbH & Co. KG