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# Processing instructions according to EN ISO 17664:2017

**Product: Dental Mirror**

**Medical device**

Durotanic No.4 [REF] 355 453 / Durotanic No.5 [REF] 355 454

**Manufactured by: SHAHAB SURGICAL INDUSTRIES**

**Shahab Road, Sialkot – Pakistan**

**Tel: +92-52-3251540, Cell: +92-301-8715510, Email: info@shahabsi.com**

**a. Application (Dental Mirror):** Shahab Surgical Dental Mirror is a Periodontics, hand held dental device; it is front surface fog free mirror. Shahab Surgical Dental Mirror is used to view the hidden side of the teeth by indirect vision. The mirror is used by the certified and qualified Dental Surgeons in the hospitals, clinics, and examination rooms. Mostly the dentists and Clinicians use it when they intend to diagnose or treatment at the backside of the teeth. The mirrors are also suitable for intraoral photographs in the oral cavity before and after treatments, as well as for controlling the parallelism. Professional use of Shahab Surgical Dental Mirror ensures a very clear image of the teeth and oral flora. In order to avoid that the mirror gets foggy, we recommend wiping with soft tissue shortly before using it.

**b. Warning:** Attention! Dental Mirror is a sensitive device, its use is strictly prohibited if dropped on the floor. After dropping the mirror can break into minor pieces even unnoticeable to human eyesight. If dropped dental mirror is used, its small pieces may be dropped in the patient mouth and may cause injuries. The product has to be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every further use.

Used Shahab Surgical dental mirror must be transported to the decontamination area for reprocessing in closed or covered containers to prevent unnecessary contamination risk Disinfection alone is not enough. Thorough cleaning and disinfection are decisive prerequisites for the effective sterilization. Metal brushes and scouring pads must not be used during manual cleaning. These materials will damage the surface and finish of the reflecting mirror. Softened tap water may be used for most rinsing however purified water should be used for final rinsing to prevent mineral deposits. Remove excess biologic soil from the reflecting mirror with a disposable wipe. Place Shahab Surgical dental mirror in a container covered with damp towels. The processing should be started as early as possible, 2 hours after use, however, at the latest, during the processing; the product should not be unnecessarily exposed to wetness or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice and/or the hospital. **Limitations in the processing:** The frequent processing has only minor effects on this product. The end of the product's useful life is mainly determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product. **Personal protective equipment:** For reasons of occupational safety and to minimize infections/cross infections, suitable personal protective equipment (protective clothes, protective gloves, protective goggles and nose/mouth protection) has to be used during the entire processing, because potential exists that the devices may carry contamination.

## c. INSTRUCTIONS:

### INITIAL TREATMENT AT THE PLACE OF USE

Remove surface contamination using a disposable cloth/paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) 2 hours after the application at the latest. The drying of residues or contamination of any kind on the product has to be avoided. When using the product Processing instructions according to EN ISO 17664:201 Product: Dental Mirror No 04, 05 Manufacturer: **Shahab Surgical Industries Pakistan** Application: The mirrors

are most suitable for intraoral photographs in the oral cavity before and after treatments, as well as for controlling the parallelism. The rhodium coated surface on one side of the Dental Mirror ensures a very clear image of the teeth and oral flora. In order to avoid that the mirror gets foggy, we recommend warming it up shortly before using it (i.e. in a water bath). 6 or 7 minutes for the first time at the place of use, do not use any aldehyde-containing or alcohol-containing agents as they may cause protein fixation.

#### **PREPARATION BEFORE CLEANING**

If possible, disassemble the product into its individual components. Visual inspection for damage and wear. Here remove the handle.

#### **CLEANING & DISINFECTION: BY MACHINERY (RECOMMENDED)**

Only perfectly maintained, checked, calibrated and released cleaning and disinfection equipment according to EN ISO 15883 with valid check mark (CE mark) may be used. Make sure that during machine cleaning and disinfection of several different products/individual components, the cleaning device is not overfilled and that the products/individual components do not contact each other. We recommend completing the cleaning and disinfection cycle according to the device manufacturer's and the cleaning agent and/or disinfectant manufacturer's specifications. Make sure that the cleaning agent is compatible with the products and disinfectants that might be used. The pH of the cleaning agent and disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Application with the following agents is recommended: Cleaning agent - Dr. Weigert, neodisher MA, Neutralizer - Dr. Weigert, neodisher Z. The instructions for use of the manufacturer of cleaning agent and disinfectant must strictly be observed. In particular, the concentrations and contact times have to be observed.

#### **Description of the recommended cleaning process:**

Pre-cleaning: The pre-cleaning should be completed in a water bath (at least drinking water quality) without other cleaning additions using a cleaning brush below the water surface to achieve both, sufficient cleaning of the products/individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; it should, however, not be less than 1 minute. Cleaning tools for pre-cleaning: Flexbrush cleaning brush, REF 605 254

Information on the cleaning and disinfection device: Type: Miele G7881

#### **Description of the recommended program:**

Cleaning program: Vario TD

Cleaning temperature: 0.2 – 56.4 °C

Hold time: 09:59 min

Disinfection temperature: 65.2 – 92.8 °C

Hold time: 09:58 min

Flushing technology: Fresh water system/demineralized water

Cleaning agent: Dr. Weigert, neodisher MA, Dosing: 30 mg (3 g/l)

Neutralizer: Dr. Weigert, neodisher Z, Dosing: 1 – 2 ml/L 8

The water should be of drinking water quality. Flushing and drying should be completed by machines. The machine cleaning should not exceed a temperature of 93 °C. It has to be ensured that visual contamination has been completely removed. If after the cleaning and disinfection process, there is still contamination visible, the cleaning & disinfection has to be repeated. Before further processing, the product must be free from any residues and dry.

#### **CLEANING: MANUAL**

The manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both, sufficient cleaning of the products/individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; it should, however, not be less than 1 minute. Afterwards, rinsing of the products/individual components with water (at least drink-

ing water quality) for at least 20 seconds.

#### **Description of the recommended cleaning process:**

Pre-cleaning: Remove surface contamination using a disposable cloth

Cleaning: Insert the instrument into the strainer bowl. Suspend the strainer bowl into the cleaning bath with cleaning solution, remove contamination by means of cleaning brushes.

Cleaning agent: ID 213, company Dürr Dental, concentration: 2 %

Cleaning tools: Flex brush cleaning brush, REF 605 254

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Drying at room temperature

The manual cleaning should not exceed a temperature of 45 °C. It must be ensured that visual contamination has been completely removed. If there is still visible contamination after the cleaning process, cleaning must be repeated.

#### **DISINFECTION: MANUAL**

For the manual disinfection, please only use approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products/individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not come into contact with each other. The instructions for use of the disinfectant manufacturer must be strictly observed. In particular, concentrations and contact times have to be observed. Use with the following agents is recommended: Dürr Dental ID 213 instrument disinfection.

#### **Description of the recommended disinfection process:**

Completion of proper cleaning preparation and manual cleaning

Disinfection: Insert the instrument into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution

Disinfectant: ID 213, company Dürr Dental, concentration: 2 %

Cleaning tools: Flex brush cleaning brush, REF 605 254

Contact time: 5 min +/- 8/9 seconds

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Drying at room temperature

Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Please observe any thermal restrictions that might apply while taking the disinfectant manufacturer's specifications into consideration. Afterwards, rinsing of the products/individual components with water (at least drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry.

#### **DRYING**

If you process the products manually, the products/individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature. If you process the products with machines, drying should be performed directly by the relevant equipment using hot air. For machine drying, a temperature of 100 °C should not be exceeded.

#### **MAINTENANCE, CONTROL AND CHECK**

The product does not make special requirements on maintenance. Reassemble disassembled products/individual components. You must always carry out a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product/individual component is not visibly clean, the entire processing has to be repeated or the product/individual component has to be properly disposed of.

#### **PACKAGING**

The product is to be packed in suitable and standardized transparent sterile packaging (sterilization bags) and sealed. Observe the instructions of the sterilization bag and sealing machine man-

ufacturers as well as the current, normative requirements. Products/individual parts that have not been sterilized in bags must be used immediately.

#### **STERILIZATION**

Only tested steam vacuum autoclaves may be used. Make sure that during the sterilization of several different products/individual components, the autoclave is not overfilled and that the products/individual components do not come into contact with each other. The following sterilization cycles may be completed: Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 15 minutes. To dry the products/individual components, the drying cycle of the autoclave should be configured accordingly.

#### **STORAGE**

To maintain the sterility, the products should be stored in standardized sterilization bags at a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed once again before their use.

#### **d. ADDITIONAL INFORMATION**

Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling. The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded. The manufacturer and the competent authority of your member state must be immediately notified of all serious incidents occurring in connection with the product. For the disposal of the products, there are regionally differing recommendations and regulations. Ask your competent disposal company for the current regulations in your region.

#### **e. CONTACT TO THE MANUFACTURER**

Tel: (+92) 52-3251540, Email: [info@shahabsi.com](mailto:info@shahabsi.com)

The instructions specified above are validated by the medical device manufacturer as suitable for the preparation of a medical device for its reuse. The processor is responsible to ensure that the actually completed processing with the equipment used, and materials personnel in the preparation facility reaches the desired result. To this end, verification and/or validation and routine monitoring of the procedure are required.

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**Distributor: HAGER & WERKEN GmbH & Co. KG**