BRAVO PORTABLE II PLUS

Handbuch / Manual / Manuel / Manuale

BRAVO PORTABLE II PLUS

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BRAVO PORTABLE II PLUS

GB Portable control box & Micro-motor - Manual

Please be sure to read this manual before operation.

1. Precaution

- Install the unit at your workplace and do not use in a dirty or very hot and humid environment. (-20°C, -40°C, 0-90% RH).
- Please take care not to drop the hand piece or drill unit. This will result in damage to ball bearing and will affect the durability of the motor.
- Please do not allow water or other liquids to spill onto and into all units.
- After replacing a bur, please make sure to close the chuck handle.
- Do not open the chuck handle while the hand piece is in operation.
- While the hand piece is not in operation, it is recommended to leave a drill bit (bur) inserted into the chuck after cleaning.
- Bravo Portable II Plus is a professional use only device, for use by professionals only.

Warning about battery usage

- The battery is a LI-ION rechargeable battery (consumable goods)
- Battery lifetime is limited and depends on working conditions.
- When it doesn't reach maximum speed, it will be time to replace batteries.
- Please dispose used batteries in accordance with the national and local laws.


2. Features

- Max speed : 30,000 RPM
- Run-time : Max 8 hrs.

	Max. speed	
SH300, SH300S, H200 M300Es	30,000 RPM	8 hrs.
MH20, MH20E	20,000 / 30,000 RPM	

- Charging time : 3 hrs.
- RPM digital display
- Rechargeable battery level display (3 steps)
- Two ways of speed control : Speed control dial & button
- Automatic power saving function
- Forward / Reverse rotation
- On / off foot switch
- Use with the standard bit (bur)


3. Part introduction

- ① Control box
- a) Power button
- b) Forward/Reverse button
- c) Hand/Foot button
- d) Speed control dial & button
- e) Adaptor connector
- f) Foot switch connector

- ② Handpiece / E-type handpiece
- ③ Handpiece stand
- ④ Control box stand (optional)
- ⑤ On/Off foot switch (optional)
- ⑥ Adaptor
- ⑦ Power cord
- ⑧ Clip (optional)
- ⑨ Handpiece holder (optional)

RPM Digital display

- ① ON/speed display (00-30)
- ② Overload warning display
- ③ Charging display
- ④ Reverse mode display
- ⑤ Foot mode display, Pause function
- ⑥ Battery level display
- ⑦ Handpiece connector

4. Installation & Operation


4-1. Battery recharging (Charging time : 3hrs.)

- 1) Connect the DC plug to the adaptor connector of control box.
- 2) While battery is charging, LED blinks. When it is fully charged, LED light remains on.
- 3) Rechargeable battery level display (3 steps)

- 1 step (low)
- 2 step (medium)
- 3 step (high)

- To prevent from overheating, please be sure to separate the adaptor from the control box. Also, please avoid operating the unit during charging hours.
- When recharging, please make sure to use the original recharge adapter and power cord.

4-2. Handpiece connection (in case of E-type micro-motor)

- 1) Insert a straight or a contra angle handpiece into the micro-motor
- 2) In case of the first use after purchasing a straight handpiece or a contra angle, please run the handpiece under no load state, then use.

4-3. Operation by Hand

- 1) Connect the handpiece to a control box.
- 2) To turn on, press and hold the power button for 1 second
- 3) Please turn the speed control knob toward clockwise from MIN to MAX.
- 4) On the digital display a green light will appear - change the direction by pushing the reverse button.

4-4. Operation by Foot

- 1) Connect a foot switch to a control box.
- 2) Turn a foot mode button will be turned on a green light on a digital display.
- 3) In the foot mode, a handpiece will be operated within the limit of set speed.

4-5. Pause Function

- 1) Holding the foot mode button over a second will pause the operation of handpiece

5. Bur(Bit) insertion and removal

- 1) To remove a bur, please turn a chuck handle counter-clockwise until it clicks to stop, and then pull out the bur.
- 2) To insert the bur, please turn the chuck handle counterclockwise until it clicks to stop.
- 3) Please be sure not to use bent, asymmetrical or damaged burs. Always use standard burs.

6. Safety protection system

- (1) Automatic power saving function: If there is no operation, the control box will be power off automatically after 3 minutes.
- (2) Overload protection function:
 - If a control box gets overloaded, and overload warning display turns on in orange light. The control box will stop after 5 seconds.

How to restart: Please set the dial to MIN to restart the unit after the overload stop.

7. Trouble shooting

Trouble	Cause	Remedy
Although a control box is fully charged, a handpiece fails to run.	Loose connection of a motor cord with the motor connector of control box.	Correctly plug the motor cord to the motor connector of control box.
	Defect in the motor cord.	Replace the motor cord.
	Worn carbon brushes.	Replace the carbon brushes.
	Trouble in the motor part.	Check the motor and repair if the motor has a problem
Motor stops suddenly while in use.	Overloaded or a chuck is open in operation.	See No.6 Safety protection system. Check the chuck and close the chuck.
Motor fails to run during a foot mode.	Foot switch fails to function	Set it on "F". Check and replace the foot switch.


In case any malfunction is found, consult with an official dealer.

8. Maintenance

8-1. Replacing carbon brushes




- 1) After turning off the power, remove the carbon housing cover and motor cord from the motor.
- 2) Take off carbon brush bolts, with a screwdriver(+) and then remove worn-out carbon brushes.
- 3) Replace with new carbon brushes carefully.
- 4) Assemble the handpiece.
- 5) After running the hand piece under no load state around 20,000 rpm for 30 minutes, you may then use.

Replacement cycle of carbon brushes will depend on working hours and loading conditions.




8-2. Replacing collet chuck

- 1) A hand cap and a spindle assembly will be disassembled by turning the hand cap counterclockwise with a chuck wrench.
 - 2) The collet chuck is disassembled by turning it counterclockwise with a chuck joint wrench and the chuck wrench.
 - 3) Replace the old collet chuck with a new one and then assemble the handpiece.
 - 4) Run the handpiece under no load state and confirm proper replacement, then use.
- Replacement cycle of collet chuck will depend on working hours and loading conditions.

8-3. Replacing ball bearings

- 1) Disassemble the hand cap and a spindle assembly from the front part of the handpiece.
- 2) If a C-ring is picked out of the front chuck housing, ball bearings and other supporting parts will be disassembled.
- 3) Assemble other supporting parts and replace to new ball bearings sequentially, install the C-ring again.
- 4) Assemble the hand cap and the spindle assembly.
- 5) Run the handpiece under no load state and confirm proper replacement, and then use.

9. Specification

Control Box	Input (Vdc)	Output (Vdc)	Power (W)	Weight (g)	Dimension		
					W	D	H
Bravo Portable II Plus	12	12	24	285	67	128	30
Handpiece	Input (V)	Torque (Ncm)	Max. Speed (rpm)	Weight (g)	Dimension		
					Ø	L	
SH300	DC12	2.4	30.000	143	22	129	
SH300S	DC12	2.4	30.000	135	23	135	
H200	DC12	2.7	30.000	138	25	119	
MH20	DC12	2.2	20.000/30.000	125	21	131	
M300Es	DC12	2.4	30.000	130	22,5	97	
MH20E	DC12	2.2	20.000	70	20	96	
Adaptor	Input			Output			
	100~240V, 50/60Hz			12Vdc / 2.0A			

10. Allowable rotation speed


This chart is a general guide for a maximum speed with conventional burs.

The allowable rotation speed can be slightly different, depending on bur materials or types.

Notice

- 1) Make sure not to extend the maximum speed specified by a manufacturer or distributor.
- 2) Large burs, even being operated below maximum speed, will vibrate. Speed should be reduced immediately to prevent damage to a bur, a handpiece or operator.
- 3) Be sure to insert the bur shank all the way, until it touches the back of the collet chuck, then tighten a collect chuck. If the bur is used not inserted completely, it is very dangerous.

Rotation speed (rpm)

11. Warranty

- We do not guarantee its quality in case of users' carelessness.

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