LARIDENT S.R.L.
Via Lamaineigra 12R (Ponte Rimassa) 16030
Tribogna (Genoa) - ITALIEN
Tel. +39 0185 939481 - Fax +39 0185 938595
info@larident.it - www.larident.it

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EN

**Manufacturer information on reprocessing
of medical devices according to DIN ISO 17664:2021**

Manufacturer: Larident S.R.L.
Via Lamaneigra 12 r
16030 Tribogna (GE), Italy

Product: Miraclean Implant Black

[REF] 355 220 - 355 224

WARNING NOTICE: WHEN USING THE SCALERS, HANDLE WITH CARE AND BE GENTLE TO AVOID TIP BREAKAGE.

Limitation of reprocessing: max 20 times.

Place of handling: clean surface with a disposable cloth or paper towel.

Storage and transport: no special requirement.

Cleaning: use cleaning agent that takes into account the instructions for cleaning the product. After cleaning, remove cleanser thoroughly under running water. When disinfecting in thermodisinfectant do not exceed 93 °C.

Disinfection: use disinfecting solutions that take into account the instructions for disinfecting the product. Remove cleanser thoroughly under running water. When disinfecting in thermodisinfectant do not exceed 93 °C.

Maintenance: no special requirement.

Packaging: the product must be packed in the package mini-grip.

Sterilisation: autoclavable vapor sterilisation at 134°C for 5 minutes or at 121°C for 15 minutes.

Control / functional check: check the possible damages and possible distortion

Storage: no special requirement.

The above-mentioned instructions have been validated as SUITABLE for preparation of medical device and its re-use by the medical device manufacturer.

It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing.

Normally, validation and routine controls are necessary.

Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.

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