**HÄGER**  
**WERKEN**

# MIRATRAY® IMPLANT

---

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Instructions for Use

Mode d'emploi

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Instruktioner

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Upute za uporabu

Käyttöohje

Návod na použitie

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Die Modellherstellung.

## **Implant impression tray with foil technique**

Miratray Implant is a universal disposable impression tray for implant impression and is used instead of an individual tray. It is suitable for dentulous and edentulous jaws. The transparent special foil at the tray bottom simplifies the selection of the tray size and the correct fit. Miratray Implant impression trays do not have to be adapted individually in the area of the abutments as the foil is simply pierced while impression taking. The foil then lies exactly against the

abutments and overflow of the impression material is avoided, which guarantees a fast, clean and precise implant impression.

### **Selection of the Implant Tray**

The adequate tray size can be selected either by conventional intraoral measurement with a measuring probe or by simple visual fitting. Sufficient space (min. 2 – 3 mm) for the impression material should be left between the impression tray and the abutments respectively the alveolar process. For dentulous jaws, the locking screws and pins have to be chosen in a way that they sit sufficiently above the occlusion surfaces of the adjacent teeth, otherwise a piercing of the foil is not or only partially possible.

### **Preparation of the Impression Tray for Impression Taking**

You should always use a tray adhesive. For a better adjustment of the tray filled with impression material in the mouth, the exact points of abutment penetration can be marked with a pen when selecting the adequate Miratray Implant tray.

### **Impression Material**

A silicone or polyether impression material with high final hardness is recommended for impression. The instructions of the respective manufacturer have to be followed when working with the impression material.

### **Disinfection of Impression**

The impression has to be disinfected after removal from the patient mouth.

### **Hazard Warning**

Use the tray with undamaged foil only. Cracks in the foil allow impression material to leak, which may result in incorrect impressions.

### **Commercial Packing**

Introduction kit with three impression trays each (S, M, L) for maxilla and mandible (dentulous and edentulous). Refill packs (6 pcs.) for the respective size.


The product is for single use only. In case of a re-use its functionality can no longer be guaranteed and will highly increase the danger of infection. The manufacturer and the competent authority of your member state must be immediately notified of all serious incidents occurring in connection with the product.

Situation for impression. Fitting of the transparent Miratray Implant impression tray. The position of the implants can be marked.

Piercing of the foil with a light, steady pressure.

The clean point of penetration of the locking pins.

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HÄGER & WERKEN GmbH & Co. KG