ALGIMIX

DE	Bedienungsanleitung
EN	Instructions for Use
FR	Mode d'emploi
ES	Instrucciones para el uso
PT	Instruções de uso
IT	Istruzioni per l'uso
SV	Bruksanvisning
DA	Brugsanvisning
FI	Käyttöohje
NO	Bruksanvisning

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
EN Instructions for Use

The REF 452 045 Alginate Mixer allows easily preparing alginate paste and similar mixtures. This machine has been designed as a simple and reliable device for the dental clinic. In order to achieve an optimal performance and a long working life, we strongly recommend you to read carefully the following operating instructions.

INSTALLATION

1. Carefully unpack your Alginate Mixer and make sure that the following accessories are inside the packing:
 - A power cord
 - A mixing bowl
 - A spatula for plaster-alginate (REF 304 317)
2. Choose a firm, stiff and well-balanced surface to place the mixer, far from heat sources, vibrations or damp.
3. Connect the power cord to a 230 V, 50/60 Hz power source.
4. Fit the bowl (B) into the machine's fixing plate (C). To do this, fit together the two parts of the fitting and slightly turn the bowl clockwise.

ALGINATE MIXER: INSTRUCTIONS FOR USE

1. Connect the machine to a power point of 230 V, 50/60 Hz. Press general switch at the rear.
2. Place the bowl on the bowl holder. Fit bowl and bowl holder firmly and adjust them vertically.
3. Turn the bowl 1/5 turn clockwise (approx 70°).
4. Make sure it is securely fitted by pulling upwards.
5. Pour alginate and water inside the bowl in the proportions indicated by the manufacturer of the alginate used.
6. Press slow speed (tortoise) to pre-beat the mix with the spatula (provided with the machine).
7. After a few seconds, pre-beating is finished. Then, press fast speed (hare) and mix the mixing pressing it with the spatula. The spatula should be placed almost flat against the wall of the bowl (Picture 1).
8. Once mixing finishes (after just a few seconds), gather the alginate mix around the wall of the bowl moving the tip of the spatula from the bottom to the upper edge of the bowl (Picture 2).
9. Alginate mixing can be now used.
10. Place it inside a tray. Remove the bowl from the machine and clean it.

IMPORTANT: The machine is fitted with two operation switches: slow (turtle) and fast (hare). To prevent any damage of the unit it is important not to press the two keys at the same time.

MAINTENANCE AND CLEANING

In order to obtain an optimal reliability and a long working life of the machine, it is strongly recommended:


- To clean carefully the mixing bowl and the spatula every time you use the machine. This will prevent dirt accumulation and possible jams of the bowl's fitting system to the motor shaft.
- To clean the outer casing of the unit always with a cloth moistened in soapy water. Do not use dissolvents or any other inflammable product.
- To use Hager & Werken original spare parts in all the maintenance operations of the Algimate Mixer.

PRECAUTIONS

- Before connecting the machine, make sure that the power source is 230 V, 50/60 Hz.
- Do not fit the bowl into the spinning plate when the mixer motor is in operation.
- Do not allow children or unskilled personnel to operate the machine.
- Choose a rigid and well-balanced surface to place the Mixer. Besides, make sure the silicon pads are correctly fitted to the working table.
- When a cycle is finished, wait until the machine motor stops completely before proceeding to remove the bowl.
- Disconnect the machine from the power source during the cleaning operations.

TECHNICAL SPECIFICATIONS

Voltage: 230 V, 50/60 Hz
Height (bowl included): 230 mm
Power: 60 W
Width: 175 mm
Depth: 230 mm
Weight: 3 kg

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HAGER & WERKEN GmbH & Co. KG · Ackerstraße 1 · 47269 Duisburg, Germany