Fit®-Kofferdam Latex Medium **REF** 355335
Fit®-Kofferdam Latex Heavy **REF** 355336

DE	2
EN	3
FR.....	4
NL.....	5
ES.....	6
PT.....	7
IT	8
SV.....	9
PL.....	10
CZ.....	11
DK	12
SL.....	13
HR	14
SK.....	15
BG	16
GR.....	17
HU	18
NO.....	19
RO.....	20

HÄGER & WERKEN GmbH & Co. KG
Ackerstraße 1
47269 Duisburg, Germany
T +49 (203) 99269-0
info@hagerwerken.de
www.hagerwerken.de

CH **REP**
Tecor-Dental AG
Ebnatstrasse 65
8200 Schaffhausen
Schweiz

DE

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EN

INSTRUCTIONS FOR USE

1. For single use only. Do not reuse.
2. Fit-Kofferdam Latex Heavy and Medium are made from natural rubber. In rare cases, natural rubber may cause symptoms of an allergic reaction, such as itching, redness, hives, swelling, fever, shortness of breath, an asthma-like reaction, low blood pressure and shock. If these symptoms occur, stop using the product immediately and seek medical advice.

Intended use:

The following instruments – Fit-Kofferdam (Non)Latex, Fit-Kofferdam Klammern, Fit-Kofferdam Zangen, Fit-Kofferdam Rahmen etc. – are designed to prevent contamination from moisture or saliva during dental treatment.

Indication:

These products are used to isolate the tooth being treated with a rubber dam and are intended for patients requiring root canal treatment or restorative dental treatment.

Contraindications:

The product is contraindicated in patients with allergies or intolerances to one or more of its ingredients. There are no age restrictions.

Intended users:

This product should be used only by qualified dental professionals.

Instructions for use:

To disinfect the Rubberdam hole template (REF 355318) before use, wipe it down with an alcohol-soaked cloth or similar. Isopropyl alcohol is recommended. The hole template is not autoclavable. All other stainless steel products should be cleaned, disinfected and sterilised in accordance with the relevant reprocessing instructions before first use and after each subsequent use.

1. Check the compatibility of the Rubberdam clip without damaging the gums. (To prevent accidental swallowing, attach some dental floss to the clip beforehand.)
2. Use the Rubberdam hole template to mark the position of the holes on the Rubberdam.
3. Select a hole size on the punch of the Rubberdam hole pliers that corresponds to the size of the affected tooth.
4. Use the Rubberdam hole punch to make holes in the rubber dam.
5. Insert the Rubberdam clip through the opening in the Rubberdam and position it parallel to the dental arch of the Rubberdam.
6. Secure the Rubberdam clip using the Rubberdam clip pliers while folding the Rubberdam.
7. Place the Rubberdam clamp on the affected tooth, taking care not to pinch the gum. Administer an anaesthetic etc. in accordance with the prescribed methods. Do not spread the product further than necessary, as this could cause damage and lead to the clips breaking.
8. Slide the Rubberdam under the Rubberdam clip and secure it to the tooth surfaces with dental floss.
9. Take care when positioning the Rubberdam frame so as not to obstruct the patient's breathing.

*There may be swabs or streaks of white powder (talcum powder) on the surface of the Rubberdam. This does not affect the quality of the Rubberdam.

ADDITIONAL INFORMATION:

All serious incidents relating to the product must be reported immediately to the manufacturer and the competent authority in your Member State. Recommendations and regulations regarding the disposal of the products vary from region to region. Please contact your local waste management authorities to find out about the current regulations in your area.

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