**Fit®-Kofferdam Non Latex Medium****REF 355395****Fit®-Kofferdam Non Latex Heavy****REF 355396**

|          |    |
|----------|----|
| DE ..... | 2  |
| EN ..... | 3  |
| FR ..... | 4  |
| NL ..... | 5  |
| ES ..... | 6  |
| PT ..... | 7  |
| IT ..... | 8  |
| SV ..... | 9  |
| PL ..... | 10 |
| CZ ..... | 11 |
| DK ..... | 12 |
| SL ..... | 13 |
| HR ..... | 14 |
| SK ..... | 15 |
| BG ..... | 16 |
| GR ..... | 17 |
| HU ..... | 18 |
| NO ..... | 19 |
| RO ..... | 20 |

HÄGER & WERKEN GmbH & Co. KG  
Ackerstraße 1  
47269 Duisburg, Germany  
T +49 (203) 99269-0  
[info@hagerwerken.de](mailto:info@hagerwerken.de)  
[www.hagerwerken.de](http://www.hagerwerken.de)

**CH REP**  
Tecor-Dental AG  
Ebnatstrasse 65  
8200 Schaffhausen  
Schweiz

## DE

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## EN

### INSTRUCTIONS FOR USE

1. For single use only. Do not reuse.
2. Fit-Kofferdam Non Latex Heavy and Medium are made of latex-free material (TPE).

**Intended use:**

The following instruments – Fit-Kofferdam (Non)Latex, Fit-Kofferdam Klammern, Fit-Kofferdam Zangen, Fit-Kofferdam Rahmen etc. – are designed to prevent contamination from moisture or saliva during dental treatment.

**Indication:**

These products are used to isolate the tooth being treated with a rubber dam and are intended for patients requiring root canal treatment or restorative dental treatment.

**Contraindications:**

The product is contraindicated in patients with allergies or intolerances to one or more of its ingredients. There are no age restrictions.

**Intended users:**

This product should be used only by qualified dental professionals.

**Instructions for use:**

To disinfect the Rubberdam perforated template (REF 355318) before use, wipe it down with an alcohol-soaked cloth or similar. Isopropyl alcohol is recommended. The perforated template is not autoclavable. All other stainless steel products should be cleaned, disinfected and sterilised in accordance with the relevant reprocessing instructions before first use and after each subsequent use.

1. Check the compatibility of the clip without damaging the gums. (To prevent accidental swallowing, attach dental floss to the clip beforehand.)
2. Use the Rubberdam hole template to mark the position of the holes on the Rubberdam.
3. Select a hole size on the punch of the Rubberdam hole punch that corresponds to the size of the affected tooth.
4. Use the Rubberdam hole punch to make holes in the rubber dam.
5. Insert the Rubberdam clip through the opening in the Rubberdam and position it parallel to the dental arch of the Rubberdam.
6. Secure the Rubberdam clip using the Rubberdam clip pliers while folding the Rubberdam.
7. Place the Rubberdam clamp on the affected tooth, taking care not to pinch the gums. Administer an anaesthetic etc. in accordance with the prescribed methods. Do not spread the product further than necessary, as this could cause damage and lead to the clamps breaking.
8. Slide the Rubberdam under the Rubberdam clamp and secure it to the tooth surfaces using dental floss.
9. Take care when positioning the Rubberdam frame so as not to obstruct the patient's breathing.

**ADDITIONAL INFORMATION:**

All serious incidents relating to the product must be reported immediately to the manufacturer and the competent authority in your Member State. Recommendations and regulations regarding the disposal of the products vary from region to region. Please contact your local waste management authorities to find out about the current regulations in your area.

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