Produkt Fit® Kofferdam Lochzange REF **355321****Produkt Fit® Kofferdam Klammerzange** REF **355333**

DE	2
EN	6
FR	10
NL	14
IT	18
ES	22
PT	26
EL	30
HR	34
BG	38
HU	42
CS	46
SL	50
SK	54
RO	58
PL	62
DK	66
SE	70
NO	74

HÄGER & WERKEN GmbH & Co. KG
Ackerstraße 1
47269 Duisburg, Germany
T +49 (203) 99269-0
info@hagerwerken.de
www.hagerwerken.de

CH REP
Tecor-Dental AG
Ebnatstrasse 65
8200 Schaffhausen
Schweiz

DE

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INSTRUCTIONS FOR USE

INTENDED USE

The following instruments – Fit-Kofferdam (Non)Latex, Fit-Kofferdam Klammern, Fit-Kofferdam Zangen, Fit-Kofferdam Rahmen etc. – are designed to prevent contamination from moisture or saliva during dental treatment.

INDICATION

These products are used to shield the tooth being treated with a rubber dam and are intended for patients who require root canal treatment or restorative dental treatment.

CONTRAINDICATIONS

Use is contraindicated in patients with allergies or intolerances to one or more ingredients. Age restrictions do not need to be taken into account.

INTENDED USER

This product should be used only by qualified dental professionals.

INSTRUCTIONS

To disinfect the Rubberdam hole template (REF 355318) before use, wipe it with an alcohol-soaked cloth or similar. Isopropanol is best. The hole template is not autoclavable. All other stainless steel products should be cleaned, disinfected and sterilised in accordance with the relevant reprocessing instructions before first use and after each subsequent use.

1. Check the compatibility of the Rubberdam clamp without damaging the gums. (To prevent accidental swallowing, attach dental floss to the clamp beforehand.)
2. Use the Rubberdam hole template to mark the position of the holes on the dam.
3. Select a hole size on the Rubberdam hole punch that corresponds to the size of the diseased tooth.
4. Punch holes in the dam using the Rubberdam hole punch.
5. Insert the Rubberdam clamp through the opening of the Rubberdam and place it parallel to the arch of the dam.
6. Secure the Rubberdam clamp using the Rubberdam clamp pliers while folding the Rubberdam.
7. Place the Rubberdam clamp on the affected tooth, taking care not to pinch the gums. Administer anaesthetic, etc., according to the prescribed methods. Do not spread the product further than necessary, as this could cause damage and lead to the clamps breaking.
8. Slide the Rubberdam under the Rubberdam clamp and secure it to the tooth surfaces with dental floss.
9. Take care when positioning the Rubberdam frame so as not to obstruct the patient's breathing.

* There may be smears or streaks of white powder (talc) on the surface of the Rubberdam. This does not affect the quality of the dam.

PROCESSING INSTRUCTIONS ACCORDING TO EN ISO 17664:2021

Warning!

The product must be cleaned, disinfected and sterilised before first use and after each subsequent use according to these reprocessing instructions. Disinfection alone is not sufficient. Thorough cleaning and disinfection are essential prerequisites for effective sterilisation. Reprocessing should be started as soon as possible, but no later than 2 hours after use. Avoid exposing the product to unnecessary moisture or humidity during the reprocessing process. Please also observe the legal regulations valid in your country as well as the hygiene regulations of the medical practice or hospital.

Limitations in reprocessing:

Frequent reprocessing has little effect on this product. The product lifetime is primarily determined by wear and damage from use. In case of doubt, the products should always be sorted out and replaced at an early stage. Decisions regarding repeated use is the sole responsibility of the user. In the event of overuse, the manufacturer does not guarantee the function, performance and safety of the product.

Personal protective equipment:

For occupational safety reasons as well as to minimise infections/cross-infections, suitable personal protective equipment (protective clothing, protective gloves, protective goggles and nose-mouth protection) must be used during the entire reprocessing procedure.

INSTRUCTIONS:

INITIAL TREATMENT AT THE POINT OF USE

Remove surface contamination with a disposable cloth/paper towel while still at the place of use. Rinse the products with water (minimum drinking water quality) no later than 2 hours after use. Do not allow residues or impurities of any kind to dry on the product. Please do not use any agents containing aldehyde or alcohol for the initial treatment at the place of use, as these may lead to protein fixation.

PREPARATION BEFORE CLEANING

Disassemble product into individual components, if possible. Visually inspect for damage and wear.

CLEANING & DISINFECTION: BY MACHINE (RECOMMENDED)

Only a properly maintained, tested, calibrated and approved cleaning and disinfection device according to EN ISO 15883 with a valid test mark (CE mark) may be used. When cleaning and disinfecting several different products/individual components by machine, ensure that the cleaning device is not overfilled and that the products/individual components do not come into contact with each other. We recommend completing the cleaning and disinfection cycles according to the specifications of the device manufacturer and the manufacturer of the cleaning agent or disinfectant. Ensure that the cleaning agent is compatible with the products and any disinfectants used. The pH of the cleaning agent and disinfectant should be between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ethers, ketones, benzines), oxidising agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogenated hydrocarbons. We recommend the following products: Cleaning agent: Dr. Weigert, neodisher MA/Neutraliser: Dr. Weigert, neodisher Z

The manufacturer's instructions for use of the cleaning agent and disinfectant must be strictly observed. In particular, note the concentrations and exposure times to be used.

Description of the recommended cleaning process:

Pre-cleaning: The pre-cleaning should be completed in a water bath (minimum drinking water quality) without other cleaning additions using a cleaning brush below the water surface to achieve both sufficient cleaning of the products/individual components while avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component but should not be less than 1 minute.

Device for pre-cleaning: Mirasuc brush, REF 605216/Flexbrush, REF 605254

Details of the cleaning and disinfecting device: Type: Miele G7881

Description of the recommended programme:

Cleaning programme: VarioTD

Cleaning temperature: 50.1 - 55.6 °C

Hold time: 11:00 min.

Disinfection temperature: 65.1 - 94.4 °C

Hold time: 10:01min.

Rinsing method: Fresh water system/demineralised water

Cleaning agent: Dr. Weigert, neodisher MA, Dosing: 30 mg (3 g/L)

Neutraliser: Dr. Weigert, neodisher Z, Dosing: 1 - 2 ml/L

Use only water that is drinking water quality. Rinsing and drying should be completed by machine. Machine cleaning should be completed at a temperature not exceeding 60 °C. Ensure that visible contamination has been completely removed. If visible contamination remains after the cleaning and disinfection process, the cleaning & disinfection must be repeated. The product must be free of any residues and dry before further reprocessing.

CLEANING: MANUAL

Manual cleaning should be carried out in a water bath (minimum drinking water quality) with the cleaning agent listed below and using a brush below the water surface in order to achieve both sufficient cleaning of the products/individual components while avoiding protein fixation and to protect the environment from contamination with spray water. The duration should depend on the level of contamination of the product/component but should not be less than 1 minute. Subsequently, rinse the products/individual components with water (minimum drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process:

Pre-cleaning: Remove surface contamination with a disposable cloth

Cleaning: Place the instrument in a strainer bowl, suspend the strainer bowl in a cleaning bath with cleaning solution, remove contamination using cleaning brushes.

Cleaning agent: ID 213, Dürr Dental, concentration: 2 %

Cleaning tools: Mirasuc brush, REF 605216/Flexbrush, REF 605254

Rinsing: Rinse with water (minimum drinking water quality) for 20 seconds

Drying: Dry at room temperature

Manual cleaning should not exceed a temperature of 45 °C. Ensure that visible contamination has been completely removed. If visible contamination remains after the cleaning process, the cleaning must be repeated.

DISINFECTION: MANUAL

For manual disinfection, please use only approved disinfectants with tested efficacy (CE mark, VAH/DGHM listed). Place the products/individual components in appropriate disinfectant baths according to the specifications of the disinfectant manufacturer. Ensure that the products are sufficiently covered and do not touch each other. The disinfectant manufacturer's instructions for use must be strictly observed. In particular, note the concentrations and exposure times to be used. We recommend the following products: Dürr Dental ID 213 Instrument disinfection, concentration 2 %

Description of the recommended disinfection process:

Carrying out proper cleaning preparation and manual cleaning

Disinfection: Place the instrument into the strainer bowl. Place the strainer bowl in the disinfection bath with disinfectant solution.

Disinfectant: ID 213, Dürr Dental, concentration: 2 %

Cleaning tool: Mirasuc brush, REF 605216/Flexbrush, REF 605254

Contact time: 5 min

Rinsing: Rinse in water (minimum drinking water quality) for 20 seconds

Drying: Dry at room temperature

Ensure that the disinfectant is compatible with the products and any cleaning agents used. The pH of the disinfectant should be between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ethers, ketones, benzines), oxidising agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogenated hydrocarbons. Observe any thermal restrictions, taking into account the specifications of the disinfectant manufacturer. Subsequently rinse the products/individual components with water (minimum drinking water quality) for at least 20 seconds. The product must be free of any residues and dry before further processing.

DRYING

If you reprocess the products manually, the products/individual components can be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature.

Drying in the machine reprocessing method should take place directly through the equipment with hot air. If machine drying, do not exceed a temperature of 100 °C.

MAINTENANCE, CONTROL AND TESTING

There are no special maintenance requirements for the product. Reassemble disassembled products/individual components. A visual inspection for contamination, damage, wear, deformation must always take place before and after the individual work steps. Damaged or corroded products must not be further used. If the product/component is not visibly clean, the entire reprocessing process must be repeated or the product/component must be disposed of properly.

PACKAGING

The product must be packed and sealed in suitable and standardised transparent sterile packaging (sterilisation pouches). Follow the instructions of the manufacturers of the sterilisation pouches and sealing machines and the current normative requirements. Products/individual components not sterilised in pouches must be used immediately.

STERILISATION

Use only tested steam-vacuum autoclaves. When sterilising several different products/individual components, ensure that the autoclave is not overfilled and that the products/individual components are not in contact with each other. The following sterilisation cycles can be completed: Steam sterilisation, 134 °C, hold time 5 minutes or steam sterilisation, 121 °C, hold time 20 minutes. To dry the products/individual components, the drying cycle of the autoclave should be set. Follow the instructions for use of the manufacturer of the autoclave.

STORAGE

To maintain sterility, store the products in standardised sterilisation pouches in a dry, clean place until use. If sterile packaging is damaged, reprocess the products again before use.

ADDITIONAL INFORMATION

Label the sterilised products according to the legal and national regulations. The recommended storage period for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and impacts during storage, transport and handling. The user must ensure within the scope of their quality management system that specified reprocessing cycles that may apply (see Limitations in reprocessing) are not exceeded. All serious incidents related to the product must be reported immediately to the manufacturer and to the competent authority of your Member State. Recommendations and regulations for the disposal of the products vary by region. Ask your competent disposal company about the current regulations in your region.

MANUFACTURER CONTACT DETAILS

Tel.: (+49) 203/99269 – 0, E-mail: info@hagerwerken.de

The above instructions have been validated by the medical device manufacturer as suitable for reprocessing a medical device for reuse. The person undertaking the reprocessing is responsible for ensuring that the reprocessing actually carried out with the equipment, materials and personnel used in the reprocessing facility achieves the desired result. Verification and/or validation and routine monitoring of the process are required for this purpose.

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