Pelotte Plus REF 254024

DE 2

EN 5

FR..... 8

NL..... 11

ES..... 14

PT 17

IT 20

SV..... 23

PL 26

CS 29

DA 32

FI 35

SK..... 38

BG 41

HU 44

NO 47

RO 50

DE

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EN

PROCESSING INSTRUCTIONS ACCORDING TO EN ISO 17664-1:2021

Warnings:

Attention! The product has to be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every further use. Disinfection alone is not enough. Thorough cleaning and disinfection are decisive prerequisites for the effective sterilization. The processing should be started as early as possible, 2 hours after use, however, at the latest. During the processing, the product should not be unnecessarily exposed to wetness or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the -medical practice and/or the hospital.

Limitations in the processing:

The frequent processing has only minor effects on this product. The end of the product's useful life is mainly -determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product.

Personal protective equipment:

For reasons of occupational safety and to minimize infections/cross infections, suitable personal protective equipment (protective clothes, protective gloves, protective goggles and nose/mouth protection) has to be used during the entire processing.

INSTRUCTIONS:

INITIAL TREATMENT AT THE PLACE OF USE

Remove surface contamination using a disposable cloth/paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) 2 hours after the application at the latest. The drying of residues or contamination of any kind on the product has to be avoided. When using the product for the first time at the place of use, do not use any aldehyde-containing or alcohol-containing agents as they may cause protein -fixation.

PREPARATION BEFORE CLEANING

If possible, disassemble the product into its individual components. Visual inspection for damage and wear. Visual inspection for damage and wear.

CLEANING & DISINFECTION: BY MACHINERY (RECOMMENDED)

Only perfectly maintained, checked, calibrated and released cleaning and disinfection equipment according to EN ISO 15883 with valid check mark (CE mark) may be used. Make sure that during machine cleaning and disinfection of several different products/individual components, the cleaning device is not overfilled and that the products/individual components do not contact each other. We recommend completing the cleaning and disinfection cycle according to the device manufacturer's and the cleaning agent and/or disinfectant manufacturer's specifications. Make sure that the cleaning agent is compatible with the products and disinfectants that might be used. The pH of the cleaning agent and disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Application with the following agents is recommended: Dr. Weigert, neodisher MA/Dr. Weigert, neodisher Z. The instructions for use of the device manufacturer and the manufacturer of cleaning agent and disinfectant must strictly be observed. Particularly the concentrations and contact times have to be observed.

Description of the recommended cleaning process:

Pre-cleaning: The pre-cleaning should be completed in a water bath (at least drinking water quality) without other cleaning additions using a cleaning brush below the water surface to achieve both, sufficient cleaning of the products/individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; it should, however, not be less than 1 minute.

Cleaning tools for pre-cleaning: Mirasuc brush, REF 605216

Information on the cleaning and disinfection device: Type: Miele G7881

Description of the recommended cleaning and disinfection:

Means: Please use the following flushing nozzles as adapter and before the cleaning and disinfection process; plug the products onto the flushing nozzles. For the exact application of the flushing nozzles, please refer to the

instructions of the relevant manufacturer. Flushing nozzles A 833, Miele & Cie KG, Gütersloh. Flushing nozzles A 816, Miele & Cie KG, Gütersloh

Cleaning program: VarioTD

Cleaning temperature: 50 - 56.43 °C

Hold time: 10:50 min.

Disinfection temperature: 65 - 94.4 °C

Hold time: 09:52 min.

Flushing technology: Fresh water system/demineralized water

Cleaning agent: Dr. Weigert, neodisher MA

Dosing: 30 g (3 g/l)

Neutralizer: Dr. Weigert, neodisher Z.

Dosing: 1 - 2 ml/l

The water should be of drinking water quality. Flushing and drying should be completed by machines. The machine cleaning should not exceed a temperature of 60 °C. It has to be ensured that visual contamination has been completely removed. If after the cleaning and disinfection process, there is still contamination visible, the cleaning and disinfection has to be repeated. Before further processing, the product must be free from any residues and dry.

CLEANING: MANUAL

The manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both, sufficient cleaning of the products/individual components avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; it should, however, not be less than 1 minute. Afterwards, rinsing of the products/individual components with water (at least drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process:

Pre-cleaning: Remove surface contamination using a disposable cloth

Cleaning: Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution, remove contamination by means of cleaning brushes

Cleaning agent: ID 213, company Dürr Dental

Concentration: 2 %

Cleaning tools: Mirasuc brush, REF 605216

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Drying at room temperature

The manual cleaning should not exceed a temperature of 45 °C. It has to be ensured that visual contamination has been completely removed. If after the cleaning process, there is still contamination visible, the cleaning has to be repeated.

DISINFECTION: MANUAL

For the manual disinfection, please only use approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products/individual components into corresponding disinfection baths according to the disinfectant manufacturer's specifications. Ensure that the products are sufficiently covered and do not contact each other. The instructions for use of the disinfectant manufacturer must be strictly observed. Particularly the concentrations and contact times have to be observed.

Application with the following agents is recommended: Dürr Dental ID 213 instrument disinfection.

Description of the recommended disinfection process:

Completion of proper cleaning preparation and manual cleaning.

Disinfection: Insert the instrument into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution.

Disinfectant: ID 213, company Dürr Dental

Concentration: 2 %

Cleaning tools: Mirasuc brush, REF 605216

Contact time: 5 min.

Rinsing: Rinse with water (at least drinking water quality) for 20 sec.

Drying: Drying at room temperature

Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Please observe thermal restrictions that might apply considering the disinfectant manufacturer's specifications.

Afterwards, rinsing of the products/individual components with water (at least drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry.

DRYING

If you process the products manually, the products/individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature. The drying in machine processing should be directly completed by the devices with hot air. With machine drying, the temperature of 100 °C should not be exceeded.

MAINTENANCE, CONTROL AND CHECK

The product does not make special requirements on maintenance. Re-assemble disassembled products/individual components. You must always carry out a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product/individual component is not visibly clean, the entire processing has to be repeated or the product/individual component has to be properly disposed of.

PACKAGING

The product is to be packed in suitable and standardized transparent sterile packaging (sterilization bags) and sealed. Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products/individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION

Only tested steam vacuum autoclaves may be used. Make sure that during the sterilization of several different products/individual components, the autoclave is not overfilled and that the products/individual components do not contact each other. The following sterilization cycles may be completed: Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 15 minutes. To dry the products/individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer's instructions for use.

STORAGE

To maintain the sterility, the products should be stored in standardized sterilization bags at a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed once again before their use.

ADDITIONAL INFORMATION

Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling. The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded. The manufacturer and the competent authority of your member state shall be immediately notified all serious incidents occurring in connection with the product. For the disposal of the products, there are regionally differing recommendations and regulations. Ask your competent disposal company for the current regulations in your region.

CONTACT TO THE MANUFACTURER

Phone: (+49) 203/99269-0, Email: info@hagerwerken.de

The instructions specified above are validated by the medical device manufacturer as suitable for the preparation of a medical device for its reuse. The processor is responsible to ensure that the actually completed processing with the equipment used, materials and personnel in the preparation facility reaches the desired result. To this end, verification and/or validation and routine monitoring of the procedure are required.

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