Matrix-Spanner Universal Nr. 1 REF **355043**
Matrix-Spanner Abgewinkelt Nr. 2 REF **355047**
Matrix-Spanner Junior (kurz) Nr. 3 REF **355048**

DE 3
EN 6
FR 9
NL 12
IT 15
ES 18
PT 21
EL 24
HR 27
BG 30
HU 33
CS 36
SL 39
SK 42
RO 45
PL 48
DA 51
SV 51
NO 57

HAGER & WERKEN GmbH & Co. KG
Ackerstraße 1
47269 Duisburg, Germany
T +49 (203) 99269-0
info@hagerwerken.de
www.hagerwerken.de

CH REP
Tecor-Dental AG
Ebnatstrasse 65
8200 Schaffhausen
Schweiz

DE

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EN

PROCESSING INSTRUCTIONS ACCORDING TO EN ISO 17664-1:2021

⚠ **Warnings: Attention!** The product must be cleaned, disinfected and sterilized according to these processing instructions before its first use as well as after every further use. Disinfection alone is not enough. Thorough cleaning and disinfection are decisive prerequisites for effective sterilization. The processing should be started as soon as possible but not later than 2 hours after use. During the processing, the product should not be unnecessarily exposed to moisture or humidity. Please also observe the legal regulations applicable in your country as well as the hygiene regulations of the medical practice and/or the hospital.

Limitations in the processing:
Frequent processing has only minor effects on this product. The end of the product's useful life is mainly determined by wear and damage caused by use. In case of doubt, the products should always be sorted out and replaced early. The user alone is responsible for making the decision regarding repeated use of the product. If the product is used too frequently, the manufacturer does not accept any warranty for the function, performance and safety of the product.

Personal protective equipment:
For reasons of occupational safety and to minimize infections/cross infections, suitable personal pro-tective equipment (protective clothes, protective gloves, protective goggles and nose/mouth protection) must be used during the entire processing.

INSTRUCTIONS:

INITIAL TREATMENT AT THE PLACE OF USE
Remove surface contamination using a disposable cloth/paper cloth when you are still at the place of use. Rinse the products with water (at least drinking water quality) no later than 2 hours after use. Do not allow residues or contamina-tion of any kind to dry on the product. When using the product for the first time at the place of use, do not use any alde-hyde-containing or alcohol-containing agents as they may cause protein fixation.

PREPARATION BEFORE CLEANING
If possible, disassemble the product into its individual components. Inspect for damage and wear.

CLEANING & DISINFECTION: BY MACHINERY (RECOMMENDED)
Only perfectly maintained, checked, calibrated and released cleaning and disinfection equipment according to EN ISO 15883 with valid check mark (CE mark) may be used. When machine cleaning and disinfecting several different products/individual components, ensure that the cleaning device is not overfilled and that the products/individual components do not contact each other. We recommend completing the cleaning and disinfection cycle according to the specifications of the manufacturer of the device and/or the cleaning agent. Make sure that the cleaning agent is compatible with the products and disinfectants that might be used. The pH of the cleaning agent and disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Application with the following agents is recommended: Dr. Weigert, neodisher MA / Dr. Weigert, neodisher Z
The instructions for use of the manufacturer of the cleaning agent and disinfectant must strictly be observed. In particular, note the concentration and contact times.

Description of the recommended cleaning process:
Pre-cleaning: The pre-cleaning should be completed in a water bath (at least drinking water quality) without other cleaning additions using a cleaning brush below the water surface to achieve both sufficient cleaning of the products/ individual components while avoiding protein fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; however, it should be no less than 1 minute.
Cleaning tool for pre-cleaning: Flexbrush cleaning brush, REF 605254
Information on the cleaning and disinfection device: Type: Miele G7881

Description of the recommended program:
Cleaning program:VarioTD
Cleaning temperature: 50.1 - 55.6 °C
Hold time: 11:00 min.
Disinfection temperature: 65.7 - 94.4 °C
Hold time: 10:01 min
Flushing technology: Fresh water system/demineralized water

Cleaning agent: Dr. Weigert, neodisher MA,
Dosing: 30 mg (3g/l)
Neutralizer: Dr. Weigert, neodisher Z,
Dosing: 1 - 2 ml/L
Aids: Insert rack for impression trays
The water should be of drinking water quality. Flushing and drying should be completed by machines. The machine cleaning should not exceed a temperature of 60 °C. Ensure that visible contamination has been completely removed. If visible contamination remains after cleaning and disinfection, repeat these steps. Before further processing, the product must be free from any residues and dry.

CLEANING: MANUAL
The manual cleaning should be completed in a water bath (at least drinking water quality) with the cleaning agent specified below using a brush below the water surface to achieve both sufficient cleaning of the products/individual components while avoiding fixation, and to protect the environment from contamination with splash water. The duration should be based on the level of contamination of the product/individual component; however, it should be no less than 1 minute. Afterwards, rinse the products/individual components with water (at least drinking water quality) for at least 20 seconds.

Description of the recommended cleaning process:
Pre-cleaning: Remove surface contamination using a disposable cloth
Cleaning: Insert the instrument into the strainer bowl. Hang the strainer bowl into the cleaning bath with cleaning solution: remove contamination by means of cleaning brushes.
Cleaning agent: ID 213, company Dürr Dental
Concentration: 2 %
Cleaning tool: Flexbrush cleaning brush, REF 605254
Rinsing: rinse with water (at least drinking water quality) for 20 sec.
Drying: dry at room temperature
The manual cleaning should not exceed a temperature of 45 °C. Ensure that visible contamination has been completely removed. If visible contamination remains after the cleaning process, cleaning must be repeated.

DISINFECTION: MANUAL
For the manual disinfection, please use only approved disinfectants with tested effectiveness (CE mark, VAH-/DGHM-listed). Put the products/individual components into corresponding disinfection baths according to the disinfectant manufacturer’s specifications. Ensure that the products are sufficiently covered and do not contact each other. The instructions for use of the disinfectant manufacturer must be strictly observed. In particular, observe the concentration and contact times. Application with the following agent is recommended: Dürr Dental ID 213 instrument disinfection

Description of the recommended disinfection process:
Completion of proper cleaning preparation and manual cleaning
Disinfection: Insert the instrument into the strainer bowl. Hang the strainer bowl into the disinfection bath with disinfection solution; remove contamination by means of cleaning brushes.
Disinfectant: ID 213, company Dürr Dental,
Concentration: 2 %
Cleaning tool: Flexbrush cleaning brush, REF 605254
Contact time: 5 min
Rinsing: rinse with water (at least drinking water quality) for 20 sec.
Drying: dry at room temperature
Make sure that the disinfectant is compatible with the products and cleaning agents that might be used. The pH of the disinfectant should lie between 5.5 and 8.5. Do not use organic solvents (e.g. alcohols, ether, ketones, petrol), oxidizing agents (e.g. peroxides), halogens (chlorine, iodine, bromine) or aromatic/halogen hydrocarbons. Please observe thermal restrictions that might apply considering the disinfectant manufacturer’s specifications.
Afterwards, rinse the products/ individual components with water (at least drinking water quality) for at least 20 seconds. Before further processing, the product must be free from any residues and dry.

DRYING
If you process the products manually, the products/individual components may be dried by blowing them off with filtered, oil-free compressed air according to DIN ISO 8573-1 (medical cleanliness class) or at room temperature. The drying in machine processing should be directly completed by the devices using hot air. With machine drying, the temperature should not exceed 100 °C.

MAINTENANCE, CONTROL AND CHECK
There are no special maintenance requirements for the product. Re-assemble disassembled products/individual components. Always perform a visual inspection for contamination, damage, wear and deformation before and after the individual work steps. Damaged or corroded products must no longer be used. If the product/individual

component is not visibly clean, the entire processing must be repeated or the product/individual component must be properly disposed of.

PACKAGING
Pack and seal the product in suitable and standardized transparent sterile packaging (sterilization bags). Observe the instructions of the sterilization bag and sealing machine manufacturers as well as the current, normative requirements. Products/individual parts that have not been sterilized in bags must be used immediately.

STERILIZATION
Only tested steam vacuum autoclaves may be used. When sterilizing several different products/individual components, ensure that the autoclave is not overfilled and that the products/individual components do not contact each other. The following sterilization cycles may be completed: Steam sterilization, 134 °C, hold time 5 minutes or steam sterilization, 121 °C, hold time 20 minutes. To dry the products/individual components, the drying cycle of the autoclave should be set. Use the autoclave manufacturer’s instructions for use.

STORAGE
To maintain the sterility, the products should be stored in standardized sterilization bags in a dry, clean place until they are used. If the sterile packaging is damaged, the products must be processed again before use.

ADDITIONAL INFORMATION
Mark the sterilized products according to the legal and national regulations. The recommended storage duration for sterile medical devices is described in standard DIN 58953-8 and depends on external influences and effects during storage, transport and handling. The user must ensure within the scope of their quality management system that specified processing cycles that might apply (see Limitations in the processing) are not exceeded. The manufacturer and the competent authority of your member state shall be immediately notified all serious incidents occurring in connection with the product. Recommendations and regulations for disposal of the products vary by region. Ask your competent disposal company for the current regulations in your region.

CONTACT TO THE MANUFACTURER
Phone: (+49) 203/99269 -0, Email: info@hagerwerken.de
The instructions specified above are validated by the medical device manufacturer as suitable for the preparation of a medical device for its reuse. The processor is responsible for ensuring that the actually completed processing with the equipment used, materials and personnel in the preparation facility achieves the desired result. To this end, verification and/or validation and routine monitoring of the procedure are required.

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